

PepPool: CEF (CD8), human

- Product code:** 3616-1
- Contents:** The 23 peptides in the pool are MHC class I-restricted T-cell epitopes from human Cytomegalovirus, Epstein Barr virus and Influenza virus (CEF) designed to stimulate T cells with a broad array of HLA types. The synthetic peptides are >95% pure.
- Applications:** The CEF peptide pool stimulates CD8+ T cells to produce e.g., IFN- γ , IL-2, and Granzyme B and is recommended as a positive control in ELISpot and FluoroSpot assays using human PBMC. The peptide pool can also be used in Flow cytometry.
- Instructions:** Sterile handling is recommended. Dissolve the lyophilized peptide pool by addition of 40 μ l DMSO to the vial. Then add 85 μ l of PBS, mix and aliquote and store at -20°C or below. This stock solution will have a concentration of 200 μ g/ml of each peptide.
- Dilute the stock solution 1:100 in cell culture medium to obtain 2 μ g/ml of each peptide in the cell culture. Use the peptide pool in ELISpot or FluoroSpot assay for stimulation of 250,000 cells per well. For use in flow cytometry, stimulate the cells over night with 2 μ g/ml peptide pool in the presence of Brefeldin A. Use the diluted peptide solution fresh.
- Storage:** Shipped at ambient temperature. Store frozen at -20°C or below upon receipt. After reconstitution, store aliquotes at -20°C or below. The peptide solution may be stored for one month at 4-8°C without effects on stability. We recommend the aliquots not be refrozen after initial use.
- Quantity:** One vial with 25 μ g of each peptide.

MABTECH

Note; for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com
Updated on 2022-11-22

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Peptides included in PepPool: CEF (CD8), human

Peptide	Sequence	Epitope Source		HLA Allele Restriction
1	VSDGGPNLY	Influenza A	PB1 (591-599)	HLA-A1
2	CTELKLSDY	Influenza A	NP (44-52)	HLA-A1
3	GLCTLVAML	EBV	BMLF1 (259-267)	HLA-A2
4	GILGFVFTL	Influenza M	Matrix 1 (58-66)	HLA-A2
5	NLVPMVATV	CMV	pp65 (495-503)	HLA-A2
6	ILRGVSAHK	Influenza A	NP (265-273)	HLA-A3
7	RVRAYTYSK	EBV	BRLF1 (148-156)	HLA-A3
8	RLRAEAQVK	EBV	EBNA3A (603-611)	HLA-A3
9	IVTDFSVIK	EBV	EBNA3B (416-424)	HLA-A11
10	ATIGTAMYK	EBV	BRLF1 (134-143)	HLA-A11
11	DYCNVLNKEF	EBV	BRLF1 (28-37)	HLA-A24
12	KTGGPIYKR	Influenza A	NP (91-99)	HLA-A68
13	TPRVTGGGAM	CMV	pp65 (415-429)	HLA-B7
14	RPPIFIRRL	EBV	EBNA3A (379-387)	HLA-B7
15	QAKWRLQTL	EBV	EBNA3A (158-166)	HLA-B8
16	FLRGRAYGL	EBV	EBNA3A (325-333)	HLA-B8
17	RAKFKQLL	EBV	BZLF1 (190-197)	HLA-B8
18	ELRSRYWAI	Influenza A	NP (380-388)	HLA-B8
19	RRIYDLIEL	EBV	EBNA3C (258-266)	HLA-B27
20	SRYWAIRTR	Influenza A	NP (383-391)	HLA-B27
21	YPLHEQHGM	EBV	EBNA3A (458-466)	HLA-B35
22	EENLLDFVRF	EBV	EBNA3C (281-290)	HLA-B44
23	EFFWDANDIY	CMV	pp65 (511-525)	HLA-B44