

Standard mix M: EYRAplex, mouse

Intended use

This datasheet is intended for use with EYRAplex panels including Standard mix M: EYRAplex, mouse. For research use only. Not for use in diagnostic procedures.

Contents

Standard mix included in EYRAplex kit	Batch*	Quantity
Standard mix M: EYRAplex, mouse	1	2 vials

**The batch number of the standard is also stated on the standard vial label.*

Two vials of lyophilized standard mix are provided, allowing the assay to be performed on separate occasions (e.g., when only half of the plate is used per run).

Standards

The table shows the concentrations of the highest standard point (S1) obtained when diluting standards according to **Preparation of standard curve**. Standard calibration is batch-specific.

Analyte	S1 Value (pg/ml)
IFN- γ	4300
IL-1 β	21000
IL-2	4800
IL-4	3300
IL-5	1400
IL-6	7000
IL-10	4000
IL-12 (p70)	4500
IL-17A	4200
TNF- α	1800

Storage

Store at 4-8 °C upon receipt.

Note:

For research use only. Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)

Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Updated on 2026-01-28

Mabtech, Inc.

USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

www.mabtech.com



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.