

Monoclonal Antibody to Human Interleukin-2

BIOTINYLATED

Antibody:	MT2A91
Product code:	3445-2-250
Quantity:	250 µg
Immunogen:	Recombinant human IL-2
Isotype:	Mouse IgG2b
Specificity:	Native and recombinant human IL-2. Native and recombinant IL-2 from cynomolgus macaque. Please visit www.mabtech.com for information on reactivity with other monkey species.
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Biotinylation:	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
Storage:	Store product at 4-8°C.
Applications:	Western blot
Protocols:	Please contact Mabtech for more information.

MABTECH

Note; for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on, 2022-11-14

www.mabtech.com