

Monoclonal Antibody to Mouse Interleukin-17A

BIOTINYLATED

Antibody:	MT2270
Product code:	3521-14-250
Quantity:	250 µg
Immunogen:	Recombinant mouse IL-17A
Isotype:	Mouse IgG1
Specificity:	Native and recombinant mouse IL-17A
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Biotinylation:	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For quantification of IL-17A in solution e.g. cell culture supernatants using ELISA and for enumeration of IL-17A producing cells using ELISpot. MT2270 is recommended as detection mAb in combination with coating mAb IL17-I (product code 3521-3).
Protocols:	Please contact Mabtech for ELISpot and ELISA protocols.

MABTECH

Note: for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on, 2022-11-14

www.mabtech.com