

Monoclonal antibody to Human Latent TGF- β 1

Antibody:	MT324
Product code:	3550-1-250
Quantity:	250 μ g
Immunogen:	Recombinant latent TGF- β 1
Isotype:	Mouse IgG1
Specificity:	Native and recombinant latent Transforming Growth Factor- β 1 (TGF- β 1). The antibody binds to the human Latency Associated Protein (LAP), which is a part of the latent TGF- β 1 complex. The antibody does not bind to human latent TGF- β 2 or - β 3.
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 μ m) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	Western blot. Recommended concentration 1 μ g/ml.
Protocols:	Please contact Mabtech for more information.

MABTECH

Note; for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com