

# Monoclonal Antibody to Human CD40

BIOTINYLATED

<b>Antibody:</b>	S2C6
<b>Product code:</b>	3601-2-1000
<b>Quantity:</b>	1000 µg
<b>Immunogen:</b>	Bladder carcinoma cell line, T24
<b>Isotype:</b>	Mouse IgG1
<b>Specificity:</b>	Human CD40
<b>Concentration:</b>	Supplied at 1 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
<b>Purification:</b>	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
<b>Biotinylation:</b>	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
<b>Storage:</b>	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
<b>Applications:</b>	Studies of antigen expression on cells and tissue sections. Studies of B cell activation, proliferation and differentiation.
<b>Protocols:</b>	For further information, please contact Mabtech.

# MABTECH

Note: for research use only.  
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

**Mabtech AB (Head Office)**  
Sweden  
Tel: +46 8 716 27 00  
mabtech@mabtech.com

**Mabtech, Inc.**  
USA  
Tel: +1 513 871-4500  
mabtech.usa@mabtech.com

Updated on, 2022-11-14

[www.mabtech.com](http://www.mabtech.com)