

# Monoclonal antibody to Human CD43

<b>Antibody:</b>	B1B6
<b>Product code:</b>	3602-1-1000
<b>Quantity:</b>	1000 µg
<b>Immunogen:</b>	Human T cell leukemia
<b>Isotype:</b>	Mouse IgG1
<b>Specificity:</b>	Human CD43
<b>Concentration:</b>	Supplied at 1 mg/ml in sterile filtered (0.2 µm) PBS. Contains no preservations.
<b>Purification:</b>	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
<b>Storage:</b>	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
<b>Applications:</b>	Studies of CD43 expression on cells and tissue sections by flow cytometry and immunostaining.
<b>Protocols:</b>	For further information, please contact Mabtech.

# MABTECH

Note; for research use only.  
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

**Mabtech AB (Head Office)**  
Sweden  
Tel: +46 8 716 27 00  
mabtech@mabtech.com

**Mabtech, Inc.**  
USA  
Tel: +1 513 871-4500  
mabtech.usa@mabtech.com