

Monoclonal antibody to Human CD28

Antibody:	CD28-A
Product code:	3608-1-50
Quantity:	50 µg
Immunogen:	Stimulated human leukocytes
Isotype:	Mouse IgG1
Specificity:	Recognises human CD28; cross-reacts with CD28 from rhesus and cynomolgus macaques.
Concentration:	Supplied at 100 µg/ml in sterile filtered (0.2 µm) phosphate buffered saline (PBS, pH 6). Contains no preservatives.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography. The endotoxin content is below 1 EU/mg.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	CD28-A can be used together with antigen to co-stimulate specific T cells. It is also possible to use with anti-CD3 polyclonal stimulation. Suitable anti-CD3 mAbs are CD3-2 for human T cells and CD3-1 for macaque T cells. The recommended working concentration of CD28-A is 0.1 µg/ml
Protocols:	For further information, please contact Mabtech.

MABTECH

Note; for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com