

# ELISA substrate: pNPP for ALP

<b>Product code:</b>	3652-P10
<b>Contents:</b>	Contains pNPP, p-Nitrophenyl Phosphate, in diethanolamine buffer.
<b>Applications:</b>	Ready-to-use substrate for ELISA assays with alkaline phosphatase.
<b>Instructions:</b>	Add the substrate after incubation with alkaline phosphatase and subsequent washing. Incubate for typically 30-90 minutes at room temperature. The substrate produces a yellow product. Read the absorbance at 405-410 nm.
<b>Storage:</b>	Shipped at ambient temperature. Store at 4-8°C upon receipt. Adjust to room temperature before use.
<b>Size:</b>	120 ml

# MABTECH

Note; for research use only.  
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

**Mabtech AB (Head Office)**  
Sweden  
Tel: +46 8 716 27 00  
mabtech@mabtech.com  
Updated on 2022-11-22

**Mabtech, Inc.**  
USA  
Tel: +1 513 871-4500  
mabtech.usa@mabtech.com



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.