

ELISpot plates, white

| | |
|----------------------|---|
| Product code: | 3654-WP-10 |
| Contents: | Contains 10 white, sterile PVDF membrane ELISpot plates (MSIPS4W10 from Millipore). |
| Applications: | ELISpot assay |
| Storage: | Shipped at ambient temperature. Store plates at room temperature upon receipt. |

MABTECH

Note; for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com
Updated on 2022-11-22

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.