

Monoclonal Antibody to Mouse IgG2a

BIOTINYLATED

Antibody:	MG2a
Product code:	3822-6-250
Quantity:	250 µg
Immunogen:	Mouse IgG
Isotype:	Rat IgG2a
Specificity:	Mouse IgG2a. No cross-reactivity (<0.1%) with mouse IgG1, IgG2b, IgG3, IgA, IgE, IgM or with IgG from rat, goat, rabbit, cow and human. Low cross-reactivity (2.5%) with IgG2c (IgG2a analogue in e.g. C57BL/6 and NOD mice).
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by ion-exchange chromatography.
Biotinylation:	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For detection of mouse IgG2a in e.g. antigen-specific ELISA.
Protocols:	Please contact Mabtech for more information.

MABTECH

Note; for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on, 2022-11-14

www.mabtech.com