

# Anti-human IgA mAb (MT20), ALP

<b>Product code:</b>	3860-9A-1000
<b>Antibody:</b>	MT20
<b>Enzyme:</b>	Highly purified enzyme from calf intestine with high activity especially suited for enzyme immunoassay purposes.
<b>Conjugation:</b>	The preparation has been thoroughly depleted of free ligand and free enzyme. The conjugate has high residual enzyme activity.
<b>Quantity:</b>	1 ml
<b>Applications:</b>	For quantification of IgA using e.g. ELISA. MAb MT20 is recommended as detection mAb in combination with coating mAb MT57 (product code 3860-3).
<b>Working dilution:</b>	Supplied in 0.1 M Tris buffer with 0.002% Kathon CG. Dilute 1/1000 as a general guideline. Dilution factor may vary with differing assay conditions.
<b>Storage:</b>	Shipped at ambient temperature. On arrival store product at 4-8°C.

# MABTECH

Note; for research use only.  
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

**Mabtech AB (Head Office)**  
Sweden  
Tel: +46 8 716 27 00  
mabtech@mabtech.com

**Mabtech, Inc.**  
USA  
Tel: +1 513 871-4500  
mabtech.usa@mabtech.com

Updated on 2022-11-22

[www.mabtech.com](http://www.mabtech.com)