

Monoclonal antibodies to Mouse Ig kappa and lambda chains

Antibody:	MT24/JC5-1 The reagent consists of a combination of two different monoclonal antibodies.
Product code:	3876-3-250
Quantity:	250 µg
Immunogen:	Mouse immunoglobulins and lambda light chain
Isotype:	Rat IgG1/IgG2b
Specificity:	Mouse Ig kappa and lambda light chains. MT24 is kappa chain specific and JC5-1 is lambda chain specific. No cross-reactivity (<0.1%) with Ig from rat, goat, rabbit, cow, sheep, horse or human.
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	ELISA and ELISpot
Protocols:	Please contact Mabtech for more information.

MABTECH

Note; for research use only.
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com