

# Monoclonal antibodies to Human IgM

<b>Antibody:</b>	MT11/12
<b>Product code:</b>	3880-3-250
<b>Quantity:</b>	250 µg
<b>Immunogen:</b>	Human IgM
<b>Isotype:</b>	Mouse IgG1/Mouse IgG2b
<b>Specificity:</b>	Human IgM. Cross-reacts with IgM from rhesus and cynomolgus macaques. Please visit <a href="http://www.mabtech.com">www.mabtech.com</a> for information on reactivity with other monkey species.
<b>Concentration:</b>	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
<b>Purification:</b>	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
<b>Storage:</b>	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
<b>Applications:</b>	For quantification of IgM in serum, plasma and cell culture supernatants using ELISA and for enumeration of IgM (total and antigen-specific) producing cells using ELISpot. In ELISA, MT11/12 are recommended as coating mAbs in combination with ALP-conjugated detection mAb MT22 (product code 3880-9A). In ELISpot, MT11/12 are recommended as coating mAbs in combination with biotinylated detection mAb MT22 (product code 3880-6).
<b>Protocols:</b>	Please contact Mabtech for ELISpot and ELISA protocols..

# MABTECH

Note; for research use only.  
Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

**Mabtech AB (Head Office)**  
Sweden  
Tel: +46 8 716 27 00  
[mabtech@mabtech.com](mailto:mabtech@mabtech.com)

**Mabtech, Inc.**  
USA  
Tel: +1 513 871-4500  
[mabtech.usa@mabtech.com](mailto:mabtech.usa@mabtech.com)