

## Anti-human CCL2 mAb (detection), biotin

Product code 5311-6-100

<b>Antibody</b>	Anti-human CCL2 mAb (detection)
<b>Analyte</b>	CCL2 (MCP-1)
<b>Application</b>	Bead-based multiplex
<b>Specificity</b>	Native and recombinant human CCL2
<b>Intended use</b>	This monoclonal antibody enables specific detection of human CCL2 in bead-based multiplex assays. For research use only. Not for use in diagnostic procedures.
<b>Purification</b>	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
<b>Conjugate</b>	Biotin
<b>Quantity</b>	100 µg
<b>Concentration</b>	0.5 mg/ml
<b>Supplied in</b>	PBS with 0.02% sodium azide. Sterile-filtered (0.2 µm).
<b>Storage</b>	Store product at 4-8°C.
<b>Recommendation</b>	<p>Anti-human CCL2 mAb (detection) is recommended in combination with Anti-human CCL2 mAb (capture) (product code 5311-3) in bead-based multiplex.</p> <p>It is advised to filter the working solution of the detection mAb using a 0.2 µm filter.</p>

## Note:

For research use only. Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

**Mabtech AB (Head Office)**

Sweden  
Tel: +46 8 716 27 00  
mabtech@mabtech.com

Updated on 2026-01-22

**Mabtech, Inc.**

USA  
Tel: +1 513 871-4500  
mabtech.usa@mabtech.com

**www.mabtech.com**

Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.