

Anti-human CCL3 mAb (capture), unconjugated

Product code 5312-3-100

Antibody	Anti-human CCL3 mAb (capture)
Analyte	CCL3 (MIP-1 α)
Application	Bead-based multiplex
Specificity	Native and recombinant human CCL3
Intended use	This monoclonal antibody enables specific detection of human CCL3 in bead-based multiplex assays. For research use only. Not for use in diagnostic procedures.
Immunogen	Recombinant human CCL3
Isotype	Mouse IgG1
Purification	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Conjugate	Unconjugated
Quantity	100 μ g
Concentration	0.5 mg/ml
Supplied in	PBS with 0.02% sodium azide. Sterile-filtered (0.2 μ m).
Storage	Store product at 4-8°C.
Recommendation	Anti-human CCL3 mAb (capture) is recommended in combination with Anti-human CCL3 mAb (detection) (product code 5312-6) in bead-based multiplex.

Note:

For research use only. Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)

Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Updated on 2026-01-22

Mabtech, Inc.

USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

www.mabtech.com

Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.