

Anti-human CCL4 mAb (capture), unconjugated

Product code 5313-3-100

Antibody	Anti-human CCL4 mAb (capture)
Analyte	CCL4 (MIP-1 β)
Application	Bead-based multiplex
Specificity	Native and recombinant human CCL4
Intended use	This monoclonal antibody enables specific detection of human CCL4 in bead-based multiplex assays. For research use only. Not for use in diagnostic procedures.
Immunogen	Recombinant human CCL4 (MIP-1 β)
Isotype	Mouse IgG1
Purification	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Conjugate	Unconjugated
Quantity	100 μ g
Concentration	0.5 mg/ml
Supplied in	PBS with 0.02% sodium azide. Sterile-filtered (0.2 μ m).
Storage	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing.
Recommendation	Anti-human CCL4 mAb (capture) is recommended in combination with Anti-human CCL4 mAb (detection) (product code 5313-6) in bead-based multiplex.

Note:

For research use only. Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)

Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Updated on 2026-01-22

Mabtech, Inc.

USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

www.mabtech.com

Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.