

Anti-human CD25 mAb (capture), unconjugated

Product code 5315-3-100

Antibody	Anti-human CD25 mAb (capture)
Analyte	CD25 (IL-2R α)
Application	Bead-based multiplex
Specificity	Native and recombinant human CD25 (IL-2 receptor α)
Intended use	This monoclonal antibody enables specific detection of human CD25 (IL-2R α) in bead-based multiplex assays. For research use only. Not for use in diagnostic procedures.
Immunogen	Recombinant human soluble CD25
Isotype	Mouse IgG1
Purification	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Conjugate	Unconjugated
Quantity	100 μ g
Concentration	0.5 mg/ml
Supplied in	PBS with 0.02% sodium azide. Sterile-filtered (0.2 μ m).
Storage	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing.
Recommendation	Anti-human CD25 mAb (capture) is recommended in combination with Anti-human CD25 mAb (detection) (product code 5315-3) in bead-based multiplex.

Note:

For research use only. Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)

Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Updated on 2026-01-22

Mabtech, Inc.

USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

www.mabtech.com

Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.