

Anti-human IFN- γ mAb (detection), biotin

Product code 5318-6-100

Antibody	Anti-human IFN- γ mAb (detection)
Analyte	IFN- γ
Application	Bead-based multiplex
Specificity	Native and recombinant human IFN- γ
Intended use	This monoclonal antibody enables specific detection of human IFN- γ in bead-based multiplex assays. For research use only. Not for use in diagnostic procedures.
Purification	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Conjugate	Biotin
Quantity	100 μ g
Concentration	0.5 mg/ml
Supplied in	PBS with 0.02% sodium azide. Sterile-filtered (0.2 μ m).
Storage	Store product at 4-8°C.
Recommendation	Anti-human IFN-γ mAb (detection) is recommended in combination with Anti-human IFN-γ mAb (capture) (product code 5318-3) in bead-based multiplex. It is advised to filter the working solution of the detection mAb using a 0.2 μm filter.

Note:

For research use only. Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)

Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Updated on 2026-01-22

Mabtech, Inc.

USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

www.mabtech.com

Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.