

Anti-human IL-2 mAb (capture), unconjugated

Product code 5321-3-100

Antibody	Anti-human IL-2 mAb (capture)
Analyte	IL-2
Application	Bead-based multiplex
Specificity	Native and recombinant human IL-2
Intended use	This monoclonal antibody enables specific detection of human IL-2 in bead-based multiplex assays. For research use only. Not for use in diagnostic procedures.
Immunogen	Recombinant human IL-2
Isotype	Rat/Mouse IgG2a/IgG1
Purification	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Conjugate	Unconjugated
Quantity	100 µg
Concentration	0.5 mg/ml
Supplied in	PBS with 0.02% sodium azide. Sterile-filtered (0.2 µm).
Storage	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing.
Recommendation	Anti-human IL-2 mAb (capture) is recommended in combination with Anti-human IL-2 mAb (detection) (product code 5321-6) in bead-based multiplex.

Note:

For research use only. Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)

Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Updated on 2026-01-22

Mabtech, Inc.

USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

www.mabtech.com

Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.