

**Anti-human
Latent TGF- β 1 mAb (capture), unconjugated**

Product code 5334-3-100

Antibody	Anti-human Latent TGF- β 1 mAb (capture)
Analyte	TGF- β 1 (Latent TGF- β 1)
Application	Bead-based multiplex
Specificity	Native and recombinant latent Transforming Growth Factor- β 1 (TGF- β 1). The antibody binds to the human Latency Associated Protein (LAP), which is a part of the latent TGF- β 1 complex.
Intended use	This monoclonal antibody enables specific detection of human latent TGF- β 1 in bead-based multiplex assays. For research use only. Not for use in diagnostic procedures.
Immunogen	Recombinant human latent TGF- β 1
Isotype	Mouse IgG1
Purification	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Conjugate	Unconjugated
Quantity	100 μ g
Concentration	0.5 mg/ml
Supplied in	PBS with 0.02% sodium azide. Sterile-filtered (0.2 μ m).
Storage	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing.
Recommendation	Anti-human latent TGF- β 1 mAb (capture) is recommended in combination with Anti-human latent TGF- β 1 mAb (detection) (product code 5334-6) in bead-based multiplex.

Note:

For research use only. Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Updated on 2026-01-22

Mabtech, Inc.USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com**www.mabtech.com**

Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.