

Anti-human TNF- α mAb (detection), biotin

Product code 5335-6-100

Antibody	Anti-human TNF- α mAb (detection)
Analyte	TNF- α
Application	Bead-based multiplex
Specificity	Native and recombinant human TNF- α
Intended use	This monoclonal antibody enables specific detection of human TNF- α in bead-based multiplex assays. For research use only. Not for use in diagnostic procedures.
Purification	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Conjugate	Biotin
Quantity	100 μ g
Concentration	0.5 mg/ml
Supplied in	PBS with 0.02% sodium azide. Sterile-filtered (0.2 μ m).
Storage	Store product at 4-8°C.
Recommendation	Anti-human TNF-α mAb (detection) is recommended in combination with Anti-human TNF-α mAb (capture) (product code 5335-3) in bead-based multiplex. It is advised to filter the working solution of the detection mAb using a 0.2 μm filter.

Note:

For research use only. Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)

Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Updated on 2026-01-22

Mabtech, Inc.

USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

www.mabtech.com

Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.